


WHITE PAPER

Identifying PD-1/PD-L1 interactions

Implications for clinical research and ICI patient stratification



Author

Doroteya Raykova, PhD
Application Scientist,
Navinci Diagnostics

Background

Immune checkpoints, which are inhibitory signaling pathways that can down-modulate the immune responses of T cells, are pivotal in peripheral tissues and for maintaining immune self-tolerance. The involvement of T cells as tools in cancer treatment has emerged as an attractive therapy option, as T cell specificity, memory, and adaptability to tumor heterogeneity makes the treatment broadly applicable to different types of cancer (1). Among the many molecularly defined checkpoint proteins (2,3), two of the most studied targets are programmed cell death 1 (PD-1) receptor, also known as CD279 (cluster of differentiation 279), and its ligand PD-L1, also known as B7-H1 or CD274 (4–7), which are the main focus of this paper.

For a T cell to become activated, its T cell receptor (TCR) needs to be introduced to an antigen. Yet, priming is not complete until the CD28 co-receptor on the T cell binds to CD80 or CD86 on the antigen-presenting cell (APC) (8) (Fig. 1a). However, priming of T cells leads not only to immune response induction, but also to the initiation of an inhibitory program that, in time, can deter that same response to prevent autoimmunity. Inhibition starts with the upregulation of CTLA-4 on the surface of the T cell. CTLA-4 can outcompete the TCR co-activator CD28, generating a negative regulatory effect (9). A subsequent inhibitory signal is then initiated by the *trans*-binding of PD-1 on the T cell to its ligand PD-L1. The inhibitory function of PD-1 depends on the phosphorylation of its cytoplasmic tail, leading to the recruitment of the phosphatase SHP-2

(Fig. 1c). SHP-2 dephosphorylates its downstream targets and downregulates other signaling pathways (10,11), eventually promoting anergy, exhaustion or apoptosis in antigen-specific T cells (4,12).

Immune response can be further fine-tuned by the recently described crosstalk between the CTLA-4/CD80 and PD-1/PD-L1 inhibitory pathways through an interaction between PD-L1 and CD80 in *cis* (13,14). As a result, both the CTLA-4 and the PD-1 axes are inhibited, while CD28 co-stimulation remains active (15).

Notably, in addition to APCs and other immune cells, many tumors also express PD-L1 on their surface, which facilitates tumor evasion from the immune system (16,17). Therefore, to promote sustained T cell response and pro-inflammatory cytokine production, and to keep tumors “visible” to the immune system, immune checkpoint proteins, including PD-1 and PD-L1, are considered suitable targets for blocking with monoclonal antibodies (Table 1).

Factors affecting immune checkpoint inhibition and clinical outcomes: progress and challenges

To date, immune checkpoint inhibitors (ICIs) have significantly increased the survival rates for patients with cancers such as metastatic melanoma (19,20), non-small-cell lung cancer (NSCLC) (21,22), and renal cell carcinoma (23).




Fig. 1. T cell activation.

a) An APC presents a peptide to the T-cell via its MHC glycoprotein, and the T cell recognizes it with the help of the TCR. Co-receptor CD28 aids priming by interacting with CD80. In time, this will promote the production of cytokines by cytotoxic T cells.

b) Early inhibition begins by CTLA-4 outcompeting CD28 for the binding of CD80. PD-1 is subsequently expressed and binds to its ligand PD-L1, starting the late co-inhibitory program.

c) The PD-1/PD-L1 axis is activated when the cytoplasmic tail of PD-1 becomes phosphorylated, which recruits SHP-2 and results in downstream inhibition. Its MHC glycoprotein, and the T cell recognizes it with the help of the TCR. Co-receptor CD28 aids priming by interacting with CD80. In time, this will promote the production of cytokines by cytotoxic T cells.

Table 1. Current FDA-approved immune checkpoint inhibitors (ICIs). Modified from (18).

Agent (US brand name)	Target IC	Approved conditions
Ipilimumab (Yervoy)	CTLA-4	Melanoma, MSI-H/dMMR colorectal cancer, renal cell carcinoma (in combination with nivolumab)
Nivolumab (Opdivo)	PD-1	MSI-H or dMMR colorectal cancer, head and neck squamous cell carcinomas, hepatocellular carcinoma, melanoma, Classic Hodgkin lymphoma, non-small-cell lung carcinoma, renal cell carcinoma, urothelial cancer, small-cell lung carcinoma
Pembrolizumab (Keytruda)	PD-1	Cervical cancer, gastric cancer, head and neck squamous cell carcinomas, hepatocellular carcinoma, Classic Hodgkin lymphoma, melanoma, Merkel cell carcinoma, MSI-H or dMMR colorectal cancer, non-small-cell lung carcinoma, diffuse large B-cell lymphoma, urothelial cancer
Cemiplimab (Libtayo)	PD-1	Cutaneous squamous cell carcinoma
Atezolizumab (Tecentriq)	PD-L1	Non-small-cell lung carcinoma, urothelial cancer
Avelumab (Bavencio)	PD-L1	Merkel cell carcinoma, urothelial cancer
Durvalumab (Imfinzi)	PD-L1	Non-small-cell lung carcinoma, urothelial cancer

Although a subset of patients receiving PD-1/PD-L1 inhibitors has an impressive response to the treatment, for many, the outcomes are less positive, potentially including significant adverse effects (24). In fact, less than half of those treated with ICIs benefit from disease stability or improvement (25).

Ultimately, the reason immunotherapy does not succeed is the failure to prime CD4 and CD8 T-cell responses against the tumor (26). This can be caused various factors including lack of tumor-specific antigens (particularly in tumors with a low tumor mutational burden (TMB)), insufficient antigen presentation, suppressive signals coming from the tumor microenvironment (TME), CD8 T-cell exhaustion, and insufficient colocalization of T cells and APCs (26–28). In addition, it has recently been shown that a subset of cancer cells can express both PD-1 and PD-L1, with similar pathway activation in the tumor as observed in immune cells. Thus, treatment of such patients with ICIs would produce paradoxical effects (29). Approaches to circumvent or overcome these obstacles include the use of bispecific antibodies (bsAbs) to aid immunotherapy, and

the discovery and usage of relevant biomarkers to help define patients that will respond to treatment.

Bispecific antibodies give new hopes for treatment, emphasizing the importance of interactions

BsAbs are artificially designed antibody-based molecules that have two target-binding sites, each recognizing a different epitope or antigen (30). They are categorized as cell-bridging bsAbs, which connect immune and malignant cells, and antigen-crosslinking bsAbs that are similar to monoclonal antibody ICIs but bind two epitopes instead of one. The antigens to be crosslinked can be single targets such as HER2/HER2 or two different proteins e.g., PD-1/PD-L1, PD-1/CTLA-4, PD-L1/CTLA-4 (31). Multiple bsAbs are currently in trial (albeit almost all are still in phase I or II), in the hope that they will outperform monoclonal antibodies that have shown poor efficiency in cold tumors and have elicited resistance after treatment (31). Currently, two bsAbs, KN046 (against PD-L1/CTLA-4) and tebotelimab (against PD-1/LAG-3), are in phase III trials for solid tumor treatment (32,33).


Understanding the complex effects of the TME

With respect to biomarkers, the TME plays a pivotal role in how effective immune checkpoint therapy will be. In responders, tumors present an abundance of neoantigens, have a high load of tumor-infiltrating lymphocytes (TILs), and increased cytokine secretion, among others. In contrast, in non-responders the levels of immunosuppressive cells are high, and the levels of NK cells and activated lymphocytes are very low (34). In addition, determining TMB along with PD-L1 expression may have predictive value for identifying responders. According to the CheckMate 227 trial, a TMB of ≥ 10 mutations per megabase combined with PD-L1 expression appear to indicate longer progression-free survival in some tumors, such as NSCLC, in response to treatment with ipilimumab and nivolumab (35). Notably, in contrast to these results, a meta-analysis of ten studies with ICIs used for patients with advanced urothelial carcinoma indicated that PD-L1 expression was associated with objective response rates but not overall survival, suggesting that further investigation is needed to identify patients that will benefit from ICI therapy (36).

Current approaches to assessing patient responses to ICI therapy

Despite indications that the expression levels of a single checkpoint are insufficiently informative, and that PD-L1 abundance alone is a poor guideline of responsiveness, the leading diagnostic tool currently used to identify patients eligible for PD-1/PD-L1 inhibition is still PD-L1 immunostaining (34). When combined with an estimation of TILs, it allows patient tissues to be classified into four groups: type I (PD-L1-positive with TILs indicative of adaptive immune resistance), type II (PD-L1-negative with no TILs, indicating immune ignorance), type III (PD-L1-positive with no TILs, indicating intrinsic induction), and type IV (PD-L1-negative with TILs, indicating the role of other suppressor(s) in promoting immune tolerance) (37,38). There are several different staining options available, with varying predictive value due to the fact that PD-L1 evaluation by immunohistochemistry (IHC) is not sufficiently standardized. On the one hand, different criteria, as well as a variety of antibodies and tumor types are being used, and on the other, there is an intrinsic heterogeneity in the expression of PD-L1 across different tumors (39). The commercially available antibodies produce staining that varies considerably in intensity and patterns between different products (40).

Further variability in results can also be introduced by



the lack of standardization concerning how samples are handled and processed. The KEYNOTE 010 and other trials confirmed that archival tissue (as opposed to fresh) can be used for PD-L1 staining. Yet, the lack of unification between platforms may be an additional reason PD-L1 positivity has had such widely varying clinical significance (40–42).

The Blueprint PD-L1 IHC Assay Comparison Project, an industrial-academic collaboration, aimed to compare the analytical and clinical performance of four primary PD-L1 companion or complementary assays used in clinical trials, with the goal of establishing a cross-platform standard for PD-L1 positivity (43). Three expert pathologists reviewed 39 NSCLC formalin-fixed paraffin-embedded (FFPE) tissue samples using the FDA and EMA-approved assays outlined in Table 2. The analysis of tumor cell staining demonstrated that Ventana’s SP263 and the two Dako assays produced comparable results, whereas SP142 tended to stain fewer cells. Concerning the immune cell stain, there was higher interobserver variability across all assays (40,43). Notably, 37% of the evaluated samples showed discordance between clinical levels of PD-L1 when the assay-specific cut-off scores were applied, indicating that the classification of these patients as PD-L1-positive would depend on the assay used (43). These observations emphasize the fact that the choice of an assay and, to some degree, the user’s individual

interpretation can lead to bias towards false positives or false negatives.

The importance of interaction detection

Clinical studies report conflicting data about the association between high PD-L1 expression in tumor tissue, be it on the surface of immune or cancer cells, and clinical outcomes. PD-1 expression might be better at predicting overall survival, but high PD-L1 positivity alone is not a clear indicator of whether PD-1 will also be highly expressed or not, and PD-1/PD-L1 pathway activation in the absence of immune cells speaks for a negative prognosis (29,46–48). To a degree, the inconsistencies in the observed role of PD-L1 could be explained by the fact that PD-L1 positivity in patient samples is not directly correlated to either interaction with PD-1, or pathway activation (i.e., phosphorylation of PD-1 and effector proteins). For example, one study reported that multiple patients who tested negative for PD-L1 expression with an FDA-approved IHC assay scored high for PD-1/PD-L1 interaction. Conversely, one patient sample with a positive PD-L1 status showed a minimal number of interactions (49). PD-L1-negative tumors can still respond to ICI, as PD-L1 expression is inducible upon activation of the interferon pathway. At the same time, regardless of PD-L1 status, tumors with low infiltration of TILs are unlikely to

Table 2. US Food and Drug Administration (FDA) and European Medical Agency (EMA) approved PD-L1 assays and example cutoffs relevant for non-small-cell lung cancer (modified from (43–45)).

TC = Tumor cell score, TPS = tumor proportion score, IC = immune cell score, 1L = first line, 2L = second line, NSCLC = non-small-cell lung cancer.

Therapeutic agent (Developer)	Antibody clone used in companion/complementary assay	Interpretative scoring	IHC assay positivity (cut-off for NSCLC)	Instrument and detection system required
Nivolumab ± Ipilimumab (Bristol-Myers Squibb)	28-8 (Dako)	Tumor cell membrane	1L TC≥1%	EnVision Flex on AutostainerLink 48
Pembrolizumab (Merck)	22C3 (Dako)	Tumor cell membrane	1L TPS≥1%	EnVision Flex on AutostainerLink 48
Atezolizumab (Roche or Genentech)	SP142 (Ventana)	Tumor cell membrane Infiltrating immune cells	1L TC≥50% or IC≥10%	OptiView detection and amplification on Benchmark ULTRA
Durvalumab (AstraZeneca)	SP263 (Ventana)	Tumor cell membrane	2L TPS≥1%	OptiView detection on Benchmark ULTRA

respond to treatment (50). In line with this, when Sánchez-Magraner et al. then assessed patients with known disease outcomes, they did not observe a correlation between PD-L1 status and prognosis. Yet, they found that a low number of PD-1/PD-L1 interactions was indicative of worse overall and progression-free survival, making interaction a suitable criterion for treatment with ICI (49), and possibly with bsAbs.

In light of this, it would be useful to not only use a combination of biomarkers for predicting response and prognosis, but also to evaluate the activation status of the PD-1/PD-L1 axis (as well as other immune checkpoint pathways). A further benefit would be the addition of a stain for either immune, or tumor cells (or both), which would elucidate whether the detected interaction occurs between two immune cells, between an immune and a tumor cell, or between tumor cells. An assay that quantifies PD-1/PD-L1 interaction, for example, could therefore be of great interest to basic research, but also to the clinic for patient stratification and the approval of new companion diagnostic assays, as well as for personalized medicine, development of bsAbs and to decrease treatment-associated adverse effects. However, studies based on PD-1/PD-L1 interaction and pathway activation are almost entirely lacking. The unavailability of straightforward, sensitive and robust methods that detect interactions without disrupting tissue morphology is the most likely culprit for this.

Navinci Diagnostics' kits detect protein interactions and post-translational modifications of key players in ICI inhibition

Considering the indicated value of detecting PD-1/PD-L1 axis activation for potential patient selection and/or treatment, Navinci Diagnostics now provides a comprehensive package of products that could become an alternative to the FDA and EMA-approved PD-L1 IHC assays.

Naveni PD1/PD-L1 is a tool for direct visualization of the interaction between PD-1 and PD-L1 that has been verified in various FFPE tumor tissue sections. The proprietary Naveni technology relies on the use of a carefully selected set of Navenibodies with optimized concentrations and performance (Fig. 2). Dual target recognition is advantageous over traditional IHC, as it reduces the effects of antibody cross-reactivity and in this way improves selectivity. Proximity ligation is facilitated solely when target proteins are positioned close to one

another, ensuring sensitive detection exclusively of interacting proteins. Furthermore, the signal is augmented via rolling circle amplification (RCA), which results in strong staining even when targets are of low abundance and interactions are few. This makes the assay highly sensitive. Visualization is conducted on a regular brightfield microscope, and no additional instruments are needed.












Fig. 2. The Naveni PD1/PD-L1 assay detects the interaction between two cells expressing the receptor and its ligand.

The Navenibodies bind to their target proteins. If they are in close proximity, this allows the formation of a DNA circle amplifiable by RCA. Individual proteins are thus not able to produce signal. Detection is chromogenic and could either be HRP or AP-based, which allows flexibility in case co-staining is used. The result is visualized with standard brightfield microscopy.

Table 3 shows FFPE tissue samples from different patients with a type of NSCLC evaluated by a pathologist for the expression of relevant biomarkers. All tissues were stained with two of the approved IHC assays (SP142 and SP263, Ventana) and given a respective TC grade. Consecutive sections from the same samples were then stained with the Naveni PD1/PD-L1 kit. Notably, samples No. 101-102 show high expression of PD-L1 (based on both TC scores) but are CD8-negative, suggesting the absence of cytotoxic T cells from the tumor. Nonetheless, the high expression of CD3 is evidence for the presence of other types of lymphocytes that may express PD-1. This could explain the strong and ubiquitous interaction staining observed throughout the tissue with the Naveni kit and may be an indication that these patients would be responders to ICI therapy.

Table 3. FFPE tissue sections from 8 NSCLC patient samples were evaluated by a pathologist for various biomarkers and stained either with the standard IHC protocols for PD-L1, or with the Naveni PD1/PD-L1 kit.

No	Anatomic Site/ Pathology diagnosis	Grade Stage Type	PD-L1 (SP142)	PD-L1 (SP263)	CD70	HLA1	MPO	CD3	CD8	Naveni PD1/PD-L1 AP stain
102	Lung/Squamous cell lung carcinoma	3 IIB Malignant	TC≥ 50%	TC≥25%	+++	+++	-	+++	-	
101	Lung/Squamous cell lung carcinoma	3 IIB Malignant	TC≥ 50%	TC≥25%	+++	+++	-	+++	-	
100	Lung/Squamous cell lung carcinoma	3 IIA Malignant	TC≥25%	-	+++	-	++	-	-	
5	Lung/Squamous cell lung carcinoma	1 IB Malignant	-	-	++	+++	-	++	-	
43	Lung/Squamous cell lung carcinoma	2 IB Malignant	TC < 50% and IC < 10%	TC≥25%	+++	+++	-	+++	+	
64	Lung/Squamous cell lung carcinoma	2 IIA Malignant	TC < 50% and IC < 10%	-	+++	+++	-	+++	-	
109	Lung/Squamous cell lung carcinoma	3 IB Malignant	TC<50% and IC < 10%	TC≥25%	+++	+++	-	+++	++	
110	Lung/Squamous cell lung carcinoma	3 IB Malignant	TC<50% and IC < 10%	TC≥25%	+++	+++	-	+++	++	

(acquired from US Biomax, <https://www.biomax.us/tissue-arrays/Lung/LC1461>)
LC1461_lung_EH33-SP142_PD3-2-36

Even more interesting is sample No. 100, which is negative for lymphocyte markers (CD3-, CD8-) and antigen presentation (HLA1-) and has a moderate TC score with Ventana SP142 staining indicating some – but not high – PD-L1 expression. While sparse, the sensitive Naveni PD1/PD-L1 assay detects interactions in the tumor, suggesting the presence of both PD-1 and PD-L1 in the TME, but not on the surface of T cells. At the same time, the tissue appears to be infiltrated by neutrophils as evidenced by the relatively high expression of MPO. While neutrophils in the TME can be conditionally associated with either a favorable or an unfavorable phenotype, PD-L1+ neutrophils are known to have tumor-promoting properties as they can suppress cytotoxic T cells (51). Considering this patient’s high-grade tumor and the absence of cytotoxic cells, it is possible that the infiltrating neutrophils are expressing PD-L1. Alternatively, both PD-L1 and PD-1 may be expressed on the tumor cells (29). The detection of PD-1/PD-L1 interactions and putative pathway activation in cancer cells might indicate a patient who would not benefit from ICI treatment (29), but further studies are necessary in order to make a confident claim.

Furthermore, patient No. 109 is PD-L1-positive and presents with high immune cell infiltration, but also high expression of CD70, which is an unfavorable marker

associated with T cell exhaustion (52). There are virtually no PD-1/PD-L1 interactions detected in the tumor. In contrast, patient No. 110 presents with identical biomarker expression, but a higher activation in the PD-1/PD-L1 pathway as detected by the Naveni technology. These results suggest that a high expression of PD-L1 as indicated by IHC could be but is not necessarily evidence for the activation of the inhibitory pathway. In addition, these observations may indicate different therapeutic potential for ICI treatment in these two patients. Clinical studies that include detection of PD-1/PD-L1 interactions along with standard biomarker and IHC stains, as well as follow-up on patient survival and prognosis would be extremely valuable and have the potential to uncover more precise approaches for patient stratification.

As two alternatives, or to complement the Naveni PD1/PD-L1 kit, Naveni pY PD1 and Naveni PD1/SHP-2 sensitively detect different steps in the PD-1/PD-L1 inhibitory pathway activation in FFPE tissues. Like all Naveni-assays, Naveni pY PD1 uses dual target recognition to identify tyrosine phosphorylation of PD-1, which is the first step of PD-1-mediated inhibition (Fig. 3a and c). In comparison to traditional IHC (Fig. 3b), the assay produces a much clearer signal over background. To study the subsequent recruitment of SHP-2 to the activated PD-1 receptor,




Fig. 3. Naveni pY PD1 HRP and Naveni PD1/SHP-2 HRP performance.

Naveni pY PD1 (a) sensitively and specifically detects phosphorylated Tyr residues on PD1 in the germinal center of a human tonsil. In comparison, pY PD1 IHC staining (b) gives considerably higher background.

Consecutive FFPE sections stained with Naveni pY PD1 (c) and PD1/SHP-2 (d) consistently demonstrate the activation of the PD1/PD-L1 inhibitory pathway in the germinal center.

which is crucial in immune cells but dispensable when the pathway is activated between two tumor cells (29), one can use the upcoming Naveni PD1/SHP-2 kit (Fig. 3d). The three kits can be used on consecutive tissue sections to obtain a comprehensive picture of the pathway activation (Fig. 3c and d), which may not always correspond to PD-1 or PD-L1 expression and localization.

Conclusions

PD-1/PD-L1 immune checkpoint blockade has revolutionized the field of cancer therapy, but its efficiency has been limited to a poorly defined subset of patients. This has led to the approval of ICI companion assays which rely on IHC detection of PD-L1 in tumor tissue with the hope of identifying responders to ICI therapy. However, these assays are not standardized, fail to consider important aspects of the complex TME, and have shown an unreliable correlation to patient response and survival. To improve the predictive value of PD-L1, it is likely that different or additional biomarkers should also be evaluated(53), creating standardized biomarker panels. Furthermore, patient outcomes have a stronger correlation to high levels of PD-1/PD-L1 interaction and other hallmarks of pathway activation such as PD-1 phosphorylation and recruitment of SHP-2, than to the expression of either PD-1 or PD-L1. Navinci Diagnostics has established three assays based on proximity ligation technology that allow the reliable and sensitive detection of these three biomarkers in tumor tissue. The assays can be applied in basic research to elucidate the interplay of immune checkpoint axes and downstream molecules, in pre-clinical and clinical research to compare stainings with the existing IHC assays and evaluate the potential prognostic value of interaction detection, and in pharma, aiding the development of new drugs or bsAbs.

Available Products:

- Naveni PD1/PD-L1
- Naveni pY PD1

Coming soon:

- Naveni PD1/SHP-2

In development:

- Naveni CTLA-4/CD80
- Naveni CTLA-4/CD86
- Naveni CD28/CD80
- Naveni CD28/CD86



More information

www.navinci.se
contact@navinci.se

References

- Sharma P, Allison JP. Immune Checkpoint Targeting in Cancer Therapy: Toward Combination Strategies with Curative Potential. *Cell*. 2015 Apr;161(2):205–14.
- Zou W, Chen L. Inhibitory B7-family molecules in the tumour microenvironment. *Nat Rev Immunol*. 2008 Jun;8(6):467–77.
- Pardoll DM. The blockade of immune checkpoints in cancer immunotherapy. *Nat Rev Cancer*. 2012 Apr;12(4):252–64.
- Dong H, Strome SE, Salomao DR, Tamura H, Hirano F, Flies DB, et al. Tumor-associated B7-H1 promotes T-cell apoptosis: A potential mechanism of immune evasion. *Nat Med*. 2002 Aug 1;8(8):793–800.
- Keir ME, Butte MJ, Freeman GJ, Sharpe AH. PD-1 and Its Ligands in Tolerance and Immunity. *Annu Rev Immunol*. 2008 Apr 1;26(1):677–704.
- Sundar R, Soong R, Cho BC, Brahmer JR, Soo RA. Immunotherapy in the treatment of non-small cell lung cancer. *Lung Cancer*. 2014 Aug;85(2):101–9.
- Drake CG, Lipson EJ, Brahmer JR. Breathing new life into immunotherapy: review of melanoma, lung and kidney cancer. *Nat Rev Clin Oncol*. 2014 Jan;11(1):24–37.
- Sharma P, Wagner K, Wolchok JD, Allison JP. Novel cancer immunotherapy agents with survival benefit: recent successes and next steps. *Nat Rev Cancer*. 2011 Nov;11(11):805–12.
- Walunas TL, Lenschow DJ, Bakker CY, Linsley PS, Freeman GJ, Green JM, et al. CTLA-4 can function as a negative regulator of T cell activation. *Immunity*. 1994 Aug;1(5):405–13.
- Okazaki T, Maeda A, Nishimura H, Kurosaki T, Honjo T. PD-1 immunoreceptor inhibits B cell receptor-mediated signaling by recruiting src homology 2-domain-containing tyrosine phosphatase 2 to phosphotyrosine. *Proc Natl Acad Sci*. 2001 Nov 20;98(24):13866–71.
- Chemnitz JM, Parry RV, Nichols KE, June CH, Riley JL. SHP-1 and SHP-2 Associate with Immunoreceptor Tyrosine-Based Switch Motif of Programmed Death 1 upon Primary Human T Cell Stimulation, but Only Receptor Ligation Prevents T Cell Activation. *J Immunol*. 2004 Jul 15;173(2):945–54.
- Ostrand-Rosenberg S, Horn LA, Haile ST. The programmed death-1 immune-suppressive pathway: barrier to antitumor immunity. *J Immunol Baltim Md 1950*. 2014 Oct 15;193(8):3835–41.
- Chaudhri A, Xiao Y, Klee AN, Wang X, Zhu B, Freeman GJ. PD-L1 Binds to B7-1 Only In Cis on the Same Cell Surface. *Cancer Immunol Res*. 2018 Aug;6(8):921–9.
- Nishimura CD, Pulanco MC, Cui W, Lu L, Zang X. PD-L1 and B7-1 Cis-Interaction: New Mechanisms in Immune Checkpoints and Immunotherapies. *Trends Mol Med*. 2021 Mar;27(3):207–19.
- Zhao Y, Lee CK, Lin CH, Gassen RB, Xu X, Huang Z, et al. PD-L1:CD80 Cis-Heterodimer Triggers the Co-stimulatory Receptor CD28 While Repressing the Inhibitory PD-1 and CTLA-4 Pathways. *Immunity*. 2019 Dec 17;51(6):1059–1073.e9.
- Burr ML, Sparbier CE, Chan YC, Williamson JC, Woods K, Beavis PA, et al. CMTM6 maintains the expression of PD-L1 and regulates anti-tumour immunity. *Nature*. 2017 Sep;549(7670):101–5.
- Mezzadra R, Sun C, Jae LT, Gomez-Eerland R, de Vries E, Wu W, et al. Identification of CMTM6 and CMTM4 as PD-L1 protein regulators. *Nature*. 2017 Sep;549(7670):106–10.
- Naimi A, Mohammed RN, Raji A, Chupradit S, Yumashev AV, Suksatan W, et al. Tumor immunotherapies by immune checkpoint inhibitors (ICIs); the pros and cons. *Cell Commun Signal*. 2022 Dec;20(1):44.
- Hodi FS, O'Day SJ, McDermott DF, Weber RW, Sosman JA, Haanen JB, et al. Improved Survival with Ipilimumab in Patients with Metastatic Melanoma. *N Engl J Med*. 2010 Aug 19;363(8):711–23.
- Robert C, Ribas A, Schachter J, Arance A, Grob JJ, Mortier L, et al. Pembrolizumab versus ipilimumab in advanced melanoma (KEYNOTE-006): post-hoc 5-year results from an open-label, multicentre, randomised, controlled, phase 3 study. *Lancet Oncol*. 2019 Sep;20(9):1239–51.
- Garon EB, Hellmann MD, Rizvi NA, Carcereny E, Leigh NB, Ahn MJ, et al. Five-Year Overall Survival for Patients With Advanced Non–Small-Cell Lung Cancer Treated With Pembrolizumab: Results From the Phase I KEYNOTE-001 Study. *J Clin Oncol*. 2019 Oct 1;37(28):2518–27.
- Horn L, Spigel DR, Vokes EE, Holgado E, Ready N, Steins M, et al. Nivolumab Versus Docetaxel in Previously Treated Patients With Advanced Non–Small-Cell Lung Cancer: Two-Year Outcomes From Two Randomized, Open-Label, Phase III Trials (CheckMate 017 and CheckMate 057). *J Clin Oncol*. 2017 Dec 10;35(35):3924–33.
- Motzer RJ, Tannir NM, McDermott DF, Arén Frontera O, Melichar B, Choueiri TK, et al. Nivolumab plus Ipilimumab versus Sunitinib in Advanced Renal-Cell Carcinoma. *N Engl J Med*. 2018 Apr 5;378(14):1277–90.
- Brahmer JR, Tykodi SS, Chow LQM, Hwu WJ, Topalian SL, Hwu P, et al. Safety and Activity of Anti-PD-L1 Antibody in Patients with Advanced Cancer. *N Engl J Med*. 2012 Jun 28;366(26):2455–65.
- Haura EB, Beg AA, Rix U, Antonia S. Charting Immune Signaling Proteomes En Route to New Therapeutic Strategies. *Cancer Immunol Res*. 2015 Jul 1;3(7):714–20.
- Amelio I, Melino G, Levine AJ. Bispecific antibodies come to the aid of cancer immunotherapy. *Mol Oncol*. 2021 Jul;15(7):1759–63.
- Gupta RG, Li F, Roszik J, Lizée G. Exploiting Tumor Neoantigens to Target Cancer Evolution: Current Challenges and Promising Therapeutic Approaches. *Cancer Discov*. 2021 May;11(5):1024–39.
- Messmer MN, Snyder AG, Oberst A. Comparing the effects of different cell death programs in tumor progression and immunotherapy. *Cell Death Differ*. 2019 Jan;26(1):115–29.
- Wang X, Yang X, Zhang C, Wang Y, Cheng T, Duan L, et al. Tumor cell-intrinsic PD-1 receptor is a tumor suppressor and mediates resistance to PD-1 blockade therapy. *Proc Natl Acad Sci*. 2020 Mar 24;117(12):6640–50.
- Ma J, Mo Y, Tang M, Shen J, Qi Y, Zhao W, et al. Bispecific Antibodies: From Research to Clinical Application. *Front Immunol*. 2021 May 5;12:626616.
- Wang S, Chen K, Lei Q, Ma P, Yuan AQ, Zhao Y, et al. The state of the art of bispecific antibodies for treating human malignancies. *EMBO Mol Med [Internet]*. 2021 Sep 7 [cited 2022 Aug 31];13(9). Available from: <https://onlinelibrary.wiley.com/doi/10.15252/emmm.202114291>
- Högner A, Moehler M. Immunotherapy in Gastric Cancer. *Curr Oncol*. 2022 Mar 2;29(3):1559–74.
- Liu J, Song C, Wang X, Ni M, Wang X, Chen L, et al. Abstract P5-16-04: Preliminary safety and efficacy results of KN046 (an anti-PD-L1/CTLA-4 bispecific antibody) in combination with KN026

- (a HER2-targeted bispecific antibody) in patients with metastatic HER2-positive breast cancer: A phase II trial. *Cancer Res.* 2022 Feb 15;82(4_Supplement):P5-16-04-P5-16-04.
34. Darwin P, Toor SM, Sasidharan Nair V, Elkord E. Immune checkpoint inhibitors: recent progress and potential biomarkers. *Exp Mol Med.* 2018 Dec;50(12):1–11.
 35. Krieger T, Pearson I, Bell J, Doherty J, Robbins P. Targeted literature review on use of tumor mutational burden status and programmed cell death ligand 1 expression to predict outcomes of checkpoint inhibitor treatment. *Diagn Pathol.* 2020 Dec;15(1):6.
 36. Ghate K, Amir E, Kuksis M, Hernandez-Barajas D, Rodriguez-Romo L, Booth CM, et al. PD-L1 expression and clinical outcomes in patients with advanced urothelial carcinoma treated with checkpoint inhibitors: A meta-analysis. *Cancer Treat Rev.* 2019 Jun;76:51–6.
 37. Teng MWL, Ngiew SF, Ribas A, Smyth MJ. Classifying Cancers Based on T-cell Infiltration and PD-L1. *Cancer Res.* 2015 Jun 1;75(11):2139–45.
 38. Taube JM, Anders RA, Young GD, Xu H, Sharma R, McMiller TL, et al. Colocalization of inflammatory response with B7-h1 expression in human melanocytic lesions supports an adaptive resistance mechanism of immune escape. *Sci Transl Med.* 2012 Mar 28;4(127):127ra37.
 39. Teixidó C, Karachaliou N, González-Cao M, Morales-Espinosa D, Rosell R. Assays for predicting and monitoring responses to lung cancer immunotherapy. *Cancer Biol Med.* 2015 Jun;12(2):87–95.
 40. Grigg C, Rizvi NA. PD-L1 biomarker testing for non-small cell lung cancer: truth or fiction? *J Immunother Cancer.* 2016 Dec;4(1):48.
 41. Herbst RS, Baas P, Perez-Gracia JL, Felip E, Kim DW, Han JY, et al. Archival vs new tumor samples for assessing PD-L1 expression in the KEYNOTE-010 study of pembrolizumab (pembro) vs docetaxel (doce) for previously treated advanced NSCLC. *J Clin Oncol.* 2016 May 20;34(15_suppl):3030–3030.
 42. Herbst RS, Baas P, Perez-Gracia JL, Felip E, Kim DW, Han JY, et al. Use of archival versus newly collected tumor samples for assessing PD-L1 expression and overall survival: an updated analysis of KEYNOTE-010 trial. *Ann Oncol.* 2019 Feb;30(2):281–9.
 43. Hirsch FR, McElhinny A, Stanforth D, Ranger-Moore J, Jansson M, Kulangara K, et al. PD-L1 Immunohistochemistry Assays for Lung Cancer: Results from Phase 1 of the Blueprint PD-L1 IHC Assay Comparison Project. *J Thorac Oncol.* 2017 Feb;12(2):208–22.
 44. Jöhrens K, Rüschoff J. The Challenge to the Pathologist of PD-L1 Expression in Tumor Cells of Non-Small-Cell Lung Cancer—An Overview. *Curr Oncol.* 2021 Dec 8;28(6):5227–39.
 45. Prince EA, Sanzari JK, Pandya D, Huron D, Edwards R. Analytical Concordance of PD-L1 Assays Utilizing Antibodies From FDA-Approved Diagnostics in Advanced Cancers: A Systematic Literature Review. *JCO Precis Oncol.* 2021 Jun;5:953–73.
 46. Däster S, Eppenberger-Castori S, Mele V, Schäfer HM, Schmid L, Weixler B, et al. Low Expression of Programmed Death 1 (PD-1), PD-1 Ligand 1 (PD-L1), and Low CD8+ T Lymphocyte Infiltration Identify a Subgroup of Patients With Gastric and Esophageal Adenocarcinoma With Severe Prognosis. *Front Med.* 2020;7:144.
 47. Kurebayashi Y, Ojima H, Tsujikawa H, Kubota N, Maehara J, Abe Y, et al. Landscape of immune microenvironment in hepatocellular carcinoma and its additional impact on histological and molecular classification. *Hepatology.* 2018 Sep;68(3):1025–41.
 48. Machado I, López-Guerrero JA, Scotlandi K, Picci P, Lombart-Bosch A. Immunohistochemical analysis and prognostic significance of PD-L1, PD-1, and CD8+ tumor-infiltrating lymphocytes in Ewing's sarcoma family of tumors (ESFT). *Virchows Arch.* 2018 May;472(5):815–24.
 49. Sánchez-Magraner L, Miles J, Baker CL, Applebee CJ, Lee DJ, Elsheikh S, et al. High PD-1/PD-L1 Checkpoint Interaction Informs Tumor Selection and Therapeutic Sensitivity to Anti-PD-1/PD-L1 Treatment. *Cancer Res.* 2020 Oct 1;80(19):4244–57.
 50. Fares CM, Van Allen EM, Drake CG, Allison JP, Hu-Lieskovan S. Mechanisms of Resistance to Immune Checkpoint Blockade: Why Does Checkpoint Inhibitor Immunotherapy Not Work for All Patients? *Am Soc Clin Oncol Educ Book Am Soc Clin Oncol Annu Meet.* 2019 Jan;39:147–64.
 51. Yajuk O, Baron M, Toker S, Zelter T, Fainsod-Levi T, Granot Z. The PD-L1/PD-1 Axis Blocks Neutrophil Cytotoxicity in Cancer. *Cells.* 2021 Jun 15;10(6):1510.
 52. Ortiz-Cuaran S, Swalduz A, Foy JP, Marteau S, Morel AP, Fauvet F, et al. Epithelial-to-mesenchymal transition promotes immune escape by inducing CD70 in non-small cell lung cancer. *Eur J Cancer.* 2022 Jul;169:106–22.
 53. Wang Y, Tong Z, Zhang W, Zhang W, Buzdin A, Mu X, et al. FDA-Approved and Emerging Next Generation Predictive Biomarkers for Immune Checkpoint Inhibitors in Cancer Patients. *Front Oncol.* 2021 Jun 7;11:683419.

